
## Crown of Lamb











I. Position of rack.

- 2. Loin of lamb. Loin tail to be  $1\frac{1}{2}$  times the length of the eye muscle.
- 3. Remove the rib section of the loin.
- 4. Carefully saw through the rib bones close to the backbone.









- 5. Then remove the backbone and remove the yellow ligament (paddywack).
- 6. French trim the ribs to leave 50 mm of bone exposed. Now repeat this process to give two 7-bone racks
- 7. Make an incision between each rib. Don't cut all the way through
- 8. Attach racks to each other with twine. Ensure you attach rib ends at this point.





## Crown of Lamb – continued Loin L049 9. Form the two racks into the shape of a crown and secure loin ends with twine.



