Clod Flat Muscle (Brachialis)

Code: FQ B006

Contraction of the second seco			
I. Position of the Clod Flat Muscle.	2. To see the full cutting specification of removal from the carcase refer to specifications FQ B001, FQ B002 and FQ B003.	3. Follow the seam between the LMC and the baby LMC muscle taking care not to cut into any of the muscles	4continuing towards the natural seams of the feather muscle and remove the whole section of 4 clod muscles.
5. Once these muscles are removed take off the small thin muscle	6from the top of the baby LMC.	7. Remove the clod flat muscle.	8. Clod Fat Muscle (Brachialis).



