Silverside Mini Joints (with added fat)			Code: Silverside B003
Contraction of the second seco			
I. Position of the silverside.	2. Boneless untrimmed silverside anterior view.	3. Boneless untrimmed silverside posterior view.	4. Remove the salmon cut from the rest of the silverside by cutting along the natural seam. Remove silverwall, excess fat and connective tissue.
5. Cut into mini joints.	6. Add a 5mm layer of fat if required.	7. Secure with elasticated roasting bands.	
For this product the silverside should be matured for a minimum of I4 days.			